M-TYPER REAGENT KIT

Set 1 + Set 2

IVD

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Table of contents

1	Introduction	31
1.1	Intended use	31
1.2	Principle of the method	32
2	Kit components	32
2.1	Safety information	32
2.2	Kit components	33
2.3	Additional reagents (not included in the kit)	33
2.4	Quality control	34
2.5	Disposal of laboratory waste	34
2.6	Transport conditions	35
3	Necessary equipment	35
3.1	Auxiliary equipment	35
3.2	Recommended analytical equipment	35
4	Preparation of the measurement system	36
4.1	Preparation of the liquid chromatograph	36
4.2	Autosampler preparation	36
4.3	Mass spectrometer preparation	36
4.4	Verification of the integrated FIA-MS/MS system	39
5	Sample preparation for analysis	41
5.1	Characterization of the test material	41
5.2	Sample preparation	41
5.2.1	Sputum liquefaction	41
5.2.2	Mycolic acid isolation	41
5.2.3	Optional concentration and reconstitution step	42
5.3	Storage of the prepared sample	43
6	Sample analysis	43
6.1	Autosampler parameters	43
6.2	Liquid chromatograph program	43
6.3	Data acquisition parameters	43
6.4	Analytical sequence — recommended injection order	44
6.5	Management of abnormal signals during an analytical run	44
7	Interpretation of results	44
7.1	Principle of interpretation	44
7.2	Criteria for positive and negative results	45
7.3	Mixed samples	45
7.4	Interpretive limitations	46
8	Performance characteristics	46
8.1	Diagnostic sensitivity	46
8.2	Diagnostic specificity	46
8.3	Limit of detection	47
8.4	Other analytical parameters	47
8.5	Factors affecting analytical parameters	47
9	Limitations of the kit	47
9.1	Analytical limitations	47
9.2	Clinical limitations	48
9.3	Technical limitations	48
9.4	Storage and transport conditions affecting results	48
10	Literature	49
11	Hazardous substances	50
	Section A: Manufacturer Data	52
	Section B: Reporting of Serious Incidents	52
	Section C: Explanation of Symbols	53

1 Introduction

Mycolic acids are one of the basic structural elements of the cell wall of bacteria belonging to the suborder *Corynebacterineae*. These compounds are long-chain α -hydroxy- β -alkyl fatty acids containing two hydrocarbon chains: a longer meromycolic chain (**R₁**) and a shorter α -carbon chain. *Mycobacteria* are characterized by the presence of mycolic acids with a length of 60–90 carbon atoms while containing a fully saturated α -alkyl chain (**R₂**) with a well-defined length of 22, 24 or 26 carbon atoms.


Fig. 1. Structure of α -hydroxy- β -alkyl fatty acids.

Careful structural analyses have shown that *Mycobacterium tuberculosis* and *Mycobacterium bovis* are characterized by a definite dominance of mycolic acids with an α -alkyl chain length equal to 26 carbon atoms regardless of their type. By contrast, atypical mycobacteria are characterized by a predominance of mycolic acids with shorter alpha chains.

As early as the 1980s, researchers turned their attention to mycolic acids as potential markers for detection and differentiation of mycobacteria. These compounds are an integral part of the cell wall of mycobacteria, making up 20–40% of their dry weight. In addition, they are compounds with good chemical stability and can be detected both from material containing live mycobacteria and only their cellular sheaths or biofilm elements.

Published studies confirm the utility of mycolic acid analysis by tandem mass spectrometry in the diagnosis of tuberculosis and mycobacterioses caused by atypical mycobacteria.


Fig. 2. Types of mycolic acids.

1.1 Intended use

The ***M-Typer Reagent Kit*** is an in vitro diagnostic medical device (IVD) intended for the isolation of mycolic acids from biological material — including, but not limited to, sputum — for subsequent analysis by flow injection analysis coupled to tandem mass spectrometry (*FIA-MS/MS*). This method can serve as a supporting tool in the diagnosis of diseases caused by mycobacteria of the *Mycobacterium* genus.

Intended user: Qualified laboratory personnel in diagnostic laboratories, with practical knowledge of tandem mass spectrometer operation and MS/MS data analysis software.

Important: The test result does not constitute a standalone basis for clinical diagnosis. It should be interpreted in conjunction with other laboratory test results and the patient's clinical presentation.

1.2 Principle of the method

The ***M-Typer Reagent Kit*** is used to isolate mycolic acids from the test material. Isolation is carried out by sequentially washing off free lipids, alkaline hydrolysis, neutralization, and then liquid-liquid extraction. The obtained extract may be analyzed directly by FIA-MS/MS or stored for later analysis.

FIA-MS/MS (*flow injection analysis tandem mass spectrometry*) is based on direct injection of the sample into the chromatographic system without chromatographic separation, electrospray ionization (ESI) in negative mode, and detection by multiple reaction monitoring (MRM). Each MRM transition corresponds to a unique mycolic acid structure and enables profiling of the mycolic acid composition characteristic of a given *Mycobacterium* species.

2 Kit components

2.1 Safety information

The components of the kit contain hazardous substances. Before use, read the safety data sheets (SDS/MSDS) to assess the hazards associated with the reagents. Detailed information on each reagent can be found in the safety data sheets (SDS/MSDS) available on the manufacturer's website: www.dimedical.pl/.....

▲ SPECIAL HEALTH HAZARD WARNING (CMR): The extraction reagent (MT-004) contains substances suspected of causing cancer (H351) and suspected of damaging the unborn child (H361d). This reagent must be handled with special care. Inhalation of vapors and skin contact must be strictly avoided.

When working with the reagents, the following personal protective equipment must be used:

- disposable gloves (nitrile recommended),
- laboratory coat,
- safety glasses.

All work with reagents must be carried out under a chemical fume hood. Due to the different chemical properties of the components, the kit is supplied in two separate packages: **Set 1** and **Set 2**.

A detailed list of hazard and precautionary statements (GHS/CLP) for individual reagents is provided in the Appendix at the end of these instructions.

2.2 Kit components

All components should be stored in accordance with the conditions stated below and on the label of each reagent.

Table 1A. Kit components — Set 1.

Catalog No.	Description	Active components / Composition	Storage	Stability after opening
MT-001	Wash reagent	Main component: isopropanol. Chemical class: aliphatic alcohol. Function: washing off free lipids.	2–25°C	30 days at 2–8°C
MT-002	Hydrolysis reagent	Main component: alkaline solution of potassium hydroxide in alcohol. Chemical class: strong inorganic base in an alcoholic mixture. Function: alkaline hydrolysis of the material.	2–25°C	30 days at 2–8°C
MT-004	Extraction reagent	Main component: chloroform. Chemical class: chlorinated organic solvent. Function: liquid-liquid extraction of mycolic acids.	2–25°C	30 days at 2–8°C
MT-005	Dilution reagent	Main component: alcoholic mixture with ammonium fluoride. Chemical class: alcoholic mixture with an inorganic salt additive. Function: dilution and preparation of the extract for FIA-MS/MS analysis.	2–25°C	30 days at 2–8°C

Table 1B. List of kit components — Set 2.

Catalog No.	Description	Active components / Composition	Storage	Stability after opening
MT-003	Neutralization reagent	Main component: hydrochloric acid solution. Chemical class: inorganic acid. Function: neutralization after the hydrolysis step.	2–25°C	30 days at 2–8°C

One complete *M-Typer Reagent Kit (Set 1 + Set 2)* is intended for 100 mycolic acid isolation analyses.

The neutralization reagent (**MT-003**) is supplied in a separate package (**Set 2**) due to transport safety requirements applicable to reagents with different chemical properties. Both packages together constitute the complete kit required to perform the mycolic acid isolation procedure.

Shelf life of unopened reagents:

Unopened reagents included in the *M-Typer Reagent Kit (Set 1 and Set 2)*, when stored at 2–25°C, remain stable for 12 months from the date of manufacture.

After first opening, the reagent **must be stored at 2–8°C, kept tightly closed after each use, and used within 30 days of first opening.**

Note: Do not freeze the kit or any of its components!

2.3 Additional reagents (not included in the kit)

The following additional reagents, not supplied in the kit, are required to perform *FIA-MS/MS* measurements:

- *methanol*, of at least LC-MS/MS quality (Mobile Phase

A — MPA)

- *chloroform* stabilized with ethanol, of at least LC-MS/MS quality (Mobile Phase B — MPB and wash solvent)
- *2-propanol*, of at least LC-MS/MS quality
- 1% *N-acetyl-cysteine* in 4% *sodium hydroxide*
- sterile 0.067 mol/l phosphate buffer, pH 6.8
- 5% Tween 80 solution.

In addition, the following critical-use consumables are required:

- chromatography vials with caps and uncut Teflon septa (1.5 mL / 2 mL),
- 50 mL Falcon tubes for hydrolysis and extraction, suitable for incubation at 90°C.

2.4 Quality control

Manufacturer quality control: each kit batch is tested by the manufacturer for compliance with the quality specification before release for use.

User quality control — during routine use of the kit in the laboratory, it is recommended to:

- perform blank sample analysis at least once for every 10 analyses of real samples,
- periodically analyze a reference mycolic acid solution to verify proper performance of the *FIA-MS/MS* system,
- monitor the chromatographic signal shape and the signal-to-noise ratio.

A correct analytical signal should appear within the range of 0.05–0.30 min and should be characterized by a *hump-shaped* profile, as described in section 4.4.

If a signal is detected in the blank sample or if the signal shape is disturbed:

1. check the purity of the mobile phases,
2. flush the chromatographic system,
3. repeat the analysis.

2.5 Disposal of laboratory waste

All waste, including sealed bottles of expired reagents, is classified as chemical, hazardous waste (code 18 01 06*). Unopened packaging should be discarded in its entirety; solutions should be collected in tightly closed glass containers and submitted for disposal.

Detailed information on the disposal of individual reagents can be found in Section 13 of the corresponding *SDS/MSDS* sheets.

Test material used for analysis — regardless of its characteristics — should be considered potentially infectious and disposed of in accordance with local regulations on medical waste (code 18 01 03*).

2.6 Transport conditions

The kit should be transported at 2–25°C. Short-term temperature deviations during transport should not exceed the range of 2–30°C. Due to the chemical properties of the reagents, the ***M-Typer Reagent Kit*** is transported in two separate packages: **Set 1** and **Set 2**. Protect from light and mechanical damage.

3 Necessary equipment

3.1 Auxiliary equipment

In order to use the ***M-Typer Reagent Kit*** for the isolation of mycolic acids from biological material, the following laboratory equipment is required:

Equipment	Specification
Laboratory centrifuge	Equipped with a rotor for 50 mL tubes; minimum speed 4500 rpm. A refrigerated centrifuge (4°C) is recommended; operation at room temperature is permitted.
Vortex mixer	Standard
Horizontal laboratory shaker	For 30-minute incubation at room temperature
Heating block or water bath	Capable of incubating tubes at 90°C
Automatic pipette	Range: 100–1000 µL

3.2 Recommended analytical equipment

For *FIA-MS/MS* analysis of isolated mycolic acids, the manufacturer recommends the following analytical system:

System component	Specification
HPLC analytical pumps	Binary pump set
Autosampler	Suitable for 1.5 mL / 2 mL vials with uncut Teflon septa; injection range 1–50 µL
Tandem mass spectrometer	Type: triple quadrupole (QqQ) with ESI ion source, operating in negative ionization mode, capable of <i>FIA-MS/MS</i> analysis with MRM transition monitoring. Minimum requirements: ion detection sensitivity no lower than reference platforms (<i>Sciex QTRAP® 6500</i> , <i>Shimadzu LCMS-8050</i>). Example platforms: <i>QTRAP® 6500</i> (<i>Sciex</i>), <i>LCMS-8050</i> (<i>Shimadzu</i>) or an equivalent system meeting the above criteria.
Connecting capillary	PEEK, length 50 cm, internal diameter 100 µm

4 Preparation of the measurement system

4.1 Preparation of the liquid chromatograph

This procedure is designed for medium- and high-pressure liquid chromatographs equipped with a set of two binary pumps.

The liquid chromatograph must be prepared before analysis as follows:

1. If the chromatograph has previously worked with mobile phases other than those used in the present analysis — flush (purge) the channels of both binary pumps, first with Mobile Phase A (MPA) only.
2. After flushing with pump A — connect Mobile Phase B (MPB) to channel B and perform a purge of channel B.
3. If the chromatograph has previously worked with mobile phases according to this procedure — connect MPA and MPB directly to the respective pumps and proceed to step 4.
4. Recommended compressibility settings: pump A — 1.25 GPa; pump B — 1.60 GPa (if the chromatograph software allows).

Note: The analysis is carried out under absolutely anhydrous conditions. Prior to analysis, residual water must be removed from all components of the chromatographic system. Mobile Phase A or 2-propanol (of at least *LC-MS/MS* grade) may be used for this purpose.

4.2 Autosampler preparation

This procedure is for autosamplers allowing washing with at least one wash solution and performing injections in the range of 1–50 μL from 1.5 mL / 2 mL vials with uncut Teflon septa.

1. If the autosampler has previously worked with rinsing solutions other than those used in the present analysis — rinse (purge) it first with Mobile Phase A (MPA) only.
2. After flushing — connect the Wash Solvent and perform a purge.
3. Connect the autosampler directly to the pump and to the mass spectrometer using a PEEK capillary (50 cm / 100 μm i.d.).

Note: Deviations from the recommended capillary specification may result in differences in signal shape or shifts in the appearance times of analytical signals.

4.3 Mass spectrometer preparation

This procedure is designed for triple quadrupole (QqQ) tandem mass spectrometers equipped with an electrospray (ESI) ion source. Due to the sensitivity of the method, it is recommended to use the most sensitive instruments

available.

Prior to analysis, the detection parameters of the mass spectrometer must be optimized. The recommended parameters for the *QTRAP® 6500 (Sciex)* and *LCMS-8050 (Shimadzu)* instruments are given below. For accurate optimization of ion detection parameters, it is recommended to use a commercially available mixture of mycolic acids (e.g. Mycolic acid from *Mycobacterium tuberculosis* (bovine strain) (Product Specification [37281-34-8], Arctomsci).

Table 2A. Ion source parameters — Sciex QTRAP® 6500.

Parameter	Value
Curtain gas (CUR)	20 psi
Temperature of the ion source (TEMP)	450°C
Nebulizing gas (GS1)	50 psi
Drying gas (GS2)	70 psi
Ionization	Negative (NEG)
Voltage at the electrode (IS)	-4500 V

Table 2B. Ion source parameters — Shimadzu LCMS-8050.

Parameter	Value
Interface temperature	350°C
Ion source temperature (desolvation)	600°C
DL temperature	250°C
Nebulizing gas flow	3 L/min
Heating gas flow	10 L/min
Drying gas flow	10 L/min
Ionization mode	Negative (NEG)
Q1 resolution	Low
Q3 resolution	Unit

Table 3A. Ion transmission parameters (MRM transitions) — Sciex QTRAP® 6500.

Q1 mass [m/z]	Q3 mass [m/z]	Declustering potential [V]	Collision energy [V]
885.9	339.4	-150	-60
913.9	367.4	-150	-62
914	339.4	-150	-60
915.9	367.4	-150	-68
929.9	367.4	-150	-69
942	395.4	-150	-65
942.1	367.4	-150	-64
943.9	367.4	-150	-68
957.9	367.4	-150	-69
958	339.4	-150	-73
970	367.4	-150	-73
970	395.4	-150	-67
971.9	367.4	-150	-73
985.9	367.4	-150	-73
998	395.4	-150	-68

Q1 mass [m/z]	Q3 mass [m/z]	Declustering potential [V]	Collision energy [V]
1000	395.4	-150	-75
1014	367.4	-150	-77
1024.1	367.4	-150	-70
1028.1	395.4	-150	-77
1038.1	395.4	-150	-72
1052.1	339.4	-150	-70
1052.1	367.4	-150	-70
1080.1	367.4	-150	-74
1096.1	367.4	-150	-72
1108.1	367.4	-150	-72
1122.2	367.4	-150	-74
1124.1	395.4	-150	-74
1136.1	367.4	-150	-74
1136.2	339.4	-150	-74
1136.2	395.4	-150	-74
1150.2	367.4	-150	-74
1152.2	367.4	-150	-74
1164.2	367.4	-150	-74
1164.2	395.4	-150	-74
1166.2	367.4	-150	-74
1178.2	367.4	-150	-74
1180.2	367.4	-150	-74
1180.2	395.4	-150	-74
1192.2	367.4	-150	-74
1192.2	395.4	-150	-74
1194.2	339.4	-150	-74
1194.2	367.4	-150	-74
1196.2	395.4	-150	-74
1220.2	395.4	-150	-74
1222.2	367.4	-150	-74
1222.2	395.4	-150	-74
1250.2	367.4	-150	-74
1250.2	395.4	-150	-74
1252.2	395.4	-150	-74
1264.3	395.4	-150	-74
1280.2	395.4	-150	-74
1294.3	395.4	-150	-74

Table 3B. Ion transmission parameters (MRM transitions) — Shimadzu LCMS-8050.

Q1 mass [m/z]	Q3 mass [m/z]	Dwell time [msec]	Q1 pre bias [V]	Collision energy [V]	Q3 pre bias [V]
885.9	339.4	6	30	62	20
913.9	367.4	6	30	62	20
914	339.4	6	30	62	20
915.9	367.4	6	30	62	20
929.9	367.4	6	30	62	20
942	395.4	6	30	62	20
942.1	367.4	6	30	62	20
943.9	367.4	6	30	62	20
957.9	367.4	6	30	62	20
958	339.4	6	30	62	20
970	367.4	6	30	62	20
970	395.4	6	30	62	20
971.9	367.4	6	30	62	20
985.9	367.4	6	30	62	20
998	395.4	6	30	62	20
1000	395.4	6	30	62	20
1014	367.4	6	30	62	20
1024.1	367.4	6	30	62	20
1028.1	395.4	6	30	62	20

Q1 mass [m/z]	Q3 mass [m/z]	Dwell time [msec]	Q1 pre bias [V]	Collision energy [V]	Q3 pre bias [V]
1038.1	395.4	6	30	62	20
1052.1	339.4	6	30	62	20
1052.1	367.4	6	30	62	20
1080.1	367.4	6	30	62	20
1096.1	367.4	6	30	62	20
1108.1	367.4	6	30	62	20
1122.2	367.4	6	30	62	20
1124.1	395.4	6	30	62	20
1136.1	367.4	6	30	62	20
1136.2	339.4	6	30	62	20
1136.2	395.4	6	30	62	20
1150.2	367.4	6	30	62	20
1152.2	367.4	6	30	62	20
1164.2	367.4	6	30	62	20
1164.2	395.4	6	30	62	20
1166.2	367.4	6	30	62	20
1178.2	367.4	6	30	62	20
1180.2	367.4	6	30	62	20
1180.2	395.4	6	30	62	20
1192.2	367.4	6	30	62	20
1192.2	395.4	6	30	62	20
1194.2	339.4	6	30	62	20
1194.2	367.4	6	30	62	20
1196.2	395.4	6	30	62	20
1220.2	395.4	6	30	62	20
1222.2	367.4	6	30	62	20
1222.2	395.4	6	30	62	20
1250.2	367.4	6	30	62	20
1250.2	395.4	6	30	62	20
1252.2	395.4	6	30	62	20
1264.3	395.4	6	30	62	20
1280.2	395.4	6	30	62	20
1294.3	395.4	6	30	62	20

4.4 Verification of the integrated FIA-MS/MS system

After preparing and optimizing the parameters of the FIA-MS/MS system according to the instructions in sections 4.1–4.3, the correct operation of the system must be verified by performing *FIA-MS/MS* analysis of a solution of commercially available mycolic acids diluted in **MT-004** reagent. The analysis should be carried out in the same manner as the analysis of the test sample.

The signal from the tested substances, forming a so-called *hump shape*, should appear on the chromatogram no earlier than 0.05 min. The total duration of the signal, measured from its appearance to return to baseline, should not exceed 0.30 min. The shape of the signal should be characterized by a steep rise and fall, between which the signal remains at a relatively constant level (plateau). The plateau region should be described by at least 10 measurement points.

Deviations from the expected signal shape may result from differences in the dead volumes of the *FIA-MS/MS* system. They can be corrected by injection parameters (volume and injection rate) and the length and/or diameter of the PEEK capillary connecting the autosampler to the mass

spectrometer.


Fig. 3. Example FIA-MS/MS chromatogram showing a positive result.

Note: An additional smaller signal visible on the chromatogram after the main peak has no diagnostic significance. It is a technical artifact related to the FIA-MS/MS analysis process and the characteristics of the measurement system. Only the main analytical signal should be considered when assessing signal acceptability and interpreting the result.

5 Sample preparation for analysis

5.1 Characterization of the test material

The test material for the isolation of mycolic acids can be any material potentially containing mycobacteria of the *Mycobacterium* genus, fragments of their cell walls, or biofilm elements.

Recommended types of test material:

- sputum,
- bronchoalveolar lavage (BAL),
- cell sediment from bacterial cultures on solid or liquid media.

In justified laboratory cases, other materials may also be used, provided they have been previously validated by the laboratory..

5.2 Sample preparation

5.2.1 Sputum liquefaction

Sputum must be liquefied prior to mycolic acid isolation in accordance with the following procedure:

1. Measure **2 mL** of sputum and transfer to a 50 mL Falcon tube.
2. Add an equal volume (**2 mL**) of 1% N-acetylcysteine in 4% sodium hydroxide solution.
3. Incubate at room temperature for **20 minutes** with occasional mixing.
4. After incubation, bring the sample volume to approximately **50 mL** with sterile 0.067 mol/l phosphate buffer (pH 6.8) and add **1–2 drops** of 5% Tween 80 solution.
5. Close the tube and mix the contents by inverting several times.
6. Centrifuge at **4°C** for **20 minutes** at 4500 rpm.
7. Discard the supernatant. Subject the obtained cell pellet to further processing in accordance with section 5.2.2.

5.2.2 Mycolic acid isolation

The following procedure applies to the cell pellet obtained after sputum liquefaction (section 5.2.1) or to cell pellets obtained from other test materials.

The procedure must be carried out in the same 50 mL Falcon tube in which the cell pellet was obtained after centrifugation.

Note: The liquefaction procedure using N-acetylcysteine (section 5.2.1) applies exclusively to sputum. For material from bacterial cultures and other materials from which a cell pellet has been obtained, the mycolic acid isolation procedure must be carried out directly from the cell pellet.

Step	Action
1	Add 2 mL of wash reagent (MT-001). Vortex for 30 seconds, then centrifuge at room temperature or 4°C for 20 minutes at 4500 rpm. Discard the supernatant.
2	Add 2 mL of hydrolysis reagent (MT-002) to the washed pellet. Tightly close the tube and incubate at 90°C for 60 minutes.
3	After cooling to room temperature, add 2 mL of neutralization reagent (MT-003, Set 2). Note: a precipitate will form from the solution — this is normal.
4	Add 2 mL of extraction reagent (MT-004).
5	Perform extraction for 30 min on a horizontal shaker.
6	After phase separation, transfer the entire lower chloroform layer to a 1.5 mL autosampler vial. The volume of the transferred phase is usually approximately 1.5 mL . Note: only the lower layer should be collected for further steps. Collection of material from the upper layer may adversely affect the quality of the analytical signal. At this stage, the sample may be temporarily stored prior to further preparation for analysis, in accordance with the storage conditions specified in the instructions.
7	From the autosampler vial containing the transferred lower phase, pipette 250 µl of the sample and transfer to a new chromatographic vial. Then add 250 µl of the dilution reagent (MT-005).
8	Vortex for 30 s. The sample is ready for FIA-MS/MS analysis.

5.2.3 Optional concentration and reconstitution step

For materials with an expected low mycolic acid content, particularly paucibacillary materials such as cerebrospinal fluid, urine or other extrapulmonary materials, an additional extract concentration step may be applied after completion of mycolic acid isolation.

This step involves partial or complete evaporation of the lower chloroform phase containing the isolated mycolic acids, previously transferred to the autosampler vial in accordance with step 6 of the standard procedure. Evaporation may be carried out under a nitrogen stream or under other controlled laboratory conditions ensuring safe removal of the solvent without degradation of the analytes. After solvent evaporation, a dry residue remains in the vial. Prior to **FIA-MS/MS** analysis, the residue must be reconstituted by adding to the vial:

250 µl of reagent MT-004 — extraction reagent

and

250 µl of reagent MT-005 — dilution reagent.

After adding the reagents, vortex the vial thoroughly to ensure complete reconstitution of the residue and obtain a homogeneous final sample. The sample prepared in this way is ready for **FIA-MS/MS** analysis.

The purpose of this procedure is to increase the concentration of analytes in the final sample intended for analysis, which may improve the detectability of analytical signals in materials with low biological load.

This step is optional and should be validated by the laboratory for the given type of clinical material.

5.3 Storage of the prepared sample

The prepared sample (vial with the post-extraction fraction) may be stored prior to analysis at approximately **-70°C** for a period not exceeding **one calendar year**. Before analysis, the sample must be removed from the freezer and allowed to reach room temperature.

6 Sample analysis

FIA-MS/MS sample analysis must be performed after preparation and verification of the measurement system in accordance with sections 4.1–4.4 and after sample preparation in accordance with section 5.

6.1 Autosampler parameters

Table 4. Injection parameters.

Parameter	Value
Injection volume	50 µL
Injection rate	10 µL/s
Autosampler temperature	8°C

Note: Autosampler temperature of **8°C** is a critical parameter for sample stability and correct analysis. Exceeding this value may adversely affect the quality of the analytical signal and the reliability of the result.

6.2 Liquid chromatograph program

Table 5. Gradient program.

Time [min]	Mobile Phase A [%]	Mobile Phase B [%]	Flow [mL/min]
0.00	40	60	0.25
0.39	40	60	0.25
0.40	0	100	1.00
0.90	0	100	1.00
1.00	40	60	1.00

Total analysis time per sample: **1 minute**.

Note: The flow rate change from 0.25 mL/min to 1.00 mL/min between 0.39 and 0.40 min is a critical method parameter and must not be altered.

6.3 Data acquisition parameters

The analysis must be performed in multiple reaction monitoring (MRM) mode, in accordance with the ion transmission parameters given in Tables 3A and 3B (section 4.3).

Data acquisition requires the manufacturer’s software specific to the mass spectrometer platform used (e.g.,

Analyst for Sciex systems or **LabSolutions** for Shimadzu systems). Other MRM data acquisition software is not provided for within the approved procedure.

For processing of the obtained results, chromatographic signal integration software capable of generating reports in **TXT** or **CSV** format may be used (e.g., **MultiQuant** for Sciex systems or equivalent solutions).

Interpretation of results is performed exclusively using the dedicated **M-Typer Software**.

6.4 Analytical sequence — recommended injection order

To ensure the reliability of results, the following injection sequence is recommended in the analytical run:

1. Blank — extraction reagent (**MT-004**) transferred to an appropriate chromatography vial — verification of analytical background cleanliness.
2. Reference mycolic acid solution — verification of **FIA-MS/MS** system performance (in accordance with section 4.4).
3. Test samples.
4. Blank — analyses of test samples must be interspersed with a blank at least once every 10 analyses of real samples.

6.5 Management of abnormal signals during an analytical run

If abnormalities are detected during an analytical run, the following corrective actions should be taken:

Observation	Action
Signal in blank sample	Check the purity of mobile phases → flush the system → repeat blank
Absent signal or signal below expectation	Check injection volume and rate → check MRM parameters → repeat system verification (section 4.4)
Distorted signal shape (no plateau)	Adjust injection parameters or PEEK capillary length/diameter (section 4.2)
Signal appears outside the 0.05–0.30 min window	Check mobile phase flow → verify gradient program (Table 5)

7 Interpretation of results

7.1 Principle of interpretation

The **M-Typer Reagent Kit** is used to isolate mycolic acids from biological material. Result interpretation is performed by comparing the mycolic acid profile of the test sample with the reference profile.

Each monitored MRM transition corresponds to a signal characteristic of specific mycolic acids associated with mycobacteria of the *Mycobacterium* genus. Averaged and

normalized measurement data (mycolic acid profile) of the test sample are compared with the reference profile using the dedicated **M-Typer Software**.

M-Typer Software is a separate in vitro diagnostic medical device (IVD) bearing its own CE marking. Detailed information on result interpretation, profile classification, and taxonomic differentiation is contained in the **M-Typer Software** Instructions for Use.

Differentiation of results at the appropriate taxonomic level, including the distinction of *Mycobacterium tuberculosis* from atypical mycobacteria, requires the use of **M-Typer Software** and does not constitute an independent function of the **M-Typer Reagent Kit**.

7.2 Criteria for positive and negative results

The result of FIA-MS/MS analysis is interpreted as follows:

- **Result indicating the presence of mycobacteria of the *Mycobacterium* genus**

Obtaining a hump-shaped signal described by at least 10 measurement points in the *plateau region*, with an intensity of at least 10 times the observed background noise ($S/N \geq 10$) for a given MRM transition, measured prior to signal appearance in the interval from 0.00 to 0.05 min. A positive result is confirmed by obtaining such a signal for at least 1 of the 52 monitored MRM transitions.

- **Equivocal / indeterminate result**

Obtaining a hump-shaped signal described by at least 10 measurement points in the *plateau region*, with an intensity greater than 3 times and less than 10 times the observed background noise ($3 \leq S/N < 10$) for a given MRM transition, measured prior to signal appearance in the interval from 0.00 to 0.05 min. An equivocal result is confirmed by obtaining such a signal for at least 1 of the 52 monitored MRM transitions.

- **Result indicating the absence of mycobacteria of the *Mycobacterium* genus**

Failure to obtain a signal meeting the above criteria for any of the 52 monitored MRM transitions, i.e. absence of a *hump-shaped* signal described by at least 10 measurement points in the plateau region with an intensity exceeding 3 times the observed background noise ($S/N < 3$).

Note: Final diagnostic interpretation of the obtained mycolic acid profile, including its classification and differentiation at the appropriate taxonomic level, requires the use of the dedicated **M-Typer Software**.

7.3 Mixed samples

In the case of the presence of more than one species or strain of mycobacteria of the *Mycobacterium* genus in the test sample, the obtained mycolic acid profile may reflect a mixed signal, proportional to the quantitative contribution of the individual components of the sample. This may directly

affect the final interpretation of the result and taxonomic differentiation, which should be performed using the dedicated ***M-Typer Software***.

7.4 Interpretive limitations

The result obtained using the ***M-Typer Reagent Kit***:

- **does not constitute a standalone basis for clinical diagnosis** — it should be interpreted in conjunction with other laboratory test results and the patient's clinical presentation,
- **does not replace** bacterial culture or other standard methods of mycobacterial identification,
- is dependent on the quality of the test material, the correctness of mycolic acid isolation, and the performance of the *FIA-MS/MS* system.

8 Performance characteristics

8.1 Diagnostic sensitivity

Diagnostic sensitivity was determined on the basis of a study comprising **199 clinical specimens** from the Regional Complex of Healthcare Facilities — Centre for Pulmonary Disease Treatment and Rehabilitation in Łódź. The test material included: **131 sputum samples, 45 bronchial aspirates, 15 pleural fluids, 1 fine-needle aspiration biopsy of the lung, 5 bronchial catheter fluids, and 1 urine sample.**

Based on the combined analysis of classical, genetic, and *LC-MS/MS* results, **36 positive samples** were identified, constituting the reference point for the assessment of diagnostic sensitivity. For *LC-MS/MS*, **33 positive results** were obtained, corresponding to a diagnostic sensitivity of **92%**.

8.2 Diagnostic specificity

Diagnostic specificity was determined on the basis of the analysis of:

- 94 sputum samples collected from healthy individuals
- 88 sputum samples collected from individuals with non-mycobacterial pneumonia
- 7 bacterial culture samples from organisms outside the *Mycobacterium* genus:
 - *Rhodococcus equi*
 - *Nocardia asteroides*
 - *Corynebacterium glutamicum*
 - *Dietzia maris*
 - *Gordonia bronchialis*
 - *Gordonia sputi*
 - *Tsukamurella paurometabola*

No false-positive result was found in any of the tested samples. Diagnostic specificity is **100%**.

8.3 Limit of detection

The limit of detection (**LOD**) of the method is defined as a signal meeting the criterion $S/N \geq 3$, while the limit of quantification (**LOQ**) is defined as a signal meeting the criterion $S/N \geq 6$. The analytical sensitivity of the method was determined at 5–50 CFU/mL.

For selected strains, the following example numerical values of LOD and LOQ were determined, expressed as the number of cells per sample:

Strain	LOD [cells/sample]	LOQ [cells/sample]
<i>Mycobacterium tuberculosis</i>	6	16
<i>Mycobacterium smegmatis</i>	2	10
<i>Mycobacterium abscessus</i>	4	13

8.4 Other analytical parameters

Parameter	Value
Recovery	96–104%
Method precision (LC-MS/MS)	1–3%
Repeatability of the procedure	2–5%
Reproducibility	8–16%

Note: A full quantitative characterization of analytical interferences was not the subject of a separate validation study. From analytical practice, it is known that potential interfering signals may originate from certain plasticizers or compounds extracted from rubber components, as well as from buffer residues after the sample liquefaction step. For this reason, only the equipment, consumables, and reagents specified in sections 2.3 and 3 must be used.

8.5 Factors affecting analytical parameters

The analytical parameters of the *M-Typer Reagent Kit* depend on:

- the type and quality of the clinical material
- the correctness of the mycolic acid isolation procedure
- the technical parameters and performance of the FIA-MS/MS system
- the use of approved data analysis software

Failure to follow the kit's instructions for use may adversely affect the course and results of mycolic acid isolation and the interpretation of the test results.

9 Limitations of the kit

9.1 Analytical limitations

- The kit is intended exclusively for the isolation of mycolic acids from the test material. Due to the

potential heterogeneity and diversity of samples, the manufacturer does not define recovery parameters for individual types of clinical material.

- In the case of the presence of more than one species or strain of mycobacteria of the *Mycobacterium* genus in the test sample, the obtained mycolic acid profile will reflect averaged values, proportional to the quantitative contribution of individual species. This may lead to interpretive difficulties — see section 7.3.
- Failure to follow the instructions for use of the kit may adversely affect the course and efficiency of mycolic acid isolation and the reliability of *FIA-MS/MS* analysis results.

9.2 Clinical limitations

- The test result does not constitute a standalone basis for clinical diagnosis and must be interpreted in conjunction with other laboratory test results and the patient's clinical presentation.
- The kit is not intended to be used as the sole diagnostic tool for tuberculosis or other mycobacterioses — it is a supporting tool and does not replace bacterial culture or other standard methods of mycobacterial identification.
- The *FIA-MS/MS* method using the ***M-Typer Reagent Kit*** does not allow assessment of mycobacterial drug susceptibility.

9.3 Technical limitations

- Analysis results are dependent on the technical performance of the *FIA-MS/MS* system — including mass spectrometer calibration, purity of mobile phases, and the technical condition of the autosampler and PEEK capillary.
- The use of analytical equipment other than that recommended by the manufacturer (section 3.2) may require additional parameter optimization and may affect the analytical parameters of the kit.
- Interpretation of results requires the use of the dedicated ***M-Typer Software*** (see section 7.1).

9.4 Storage and transport conditions affecting results

- Use of reagents beyond the expiry date stated on the label is not permitted and may lead to false or unreliable results.
- Failure to observe the recommended storage or transport conditions (section 2.6) may cause degradation of reagents and a reduction in isolation performance.
- After first opening, reagents must be used within the period stated in section 2.2. Exceeding this period may affect the quality of isolation and the reliability of results.

10 Literature

The following publications constitute the scientific basis for the method and were used in the development of the M-Typer Reagent Kit.

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11 Hazardous substances

When working with reagents included in the M-Typer Reagent Kit, particular attention must be paid to the hazards listed below and appropriate precautionary measures must be taken. Detailed information on each reagent can be found in the safety data sheets (SDS/MSDS) available on the manufacturer's website: www.dimedical.pl/.....

MT-001 — Wash Reagent

Parameter	Hazard statements	
 	H225	Highly flammable liquid and vapor.
	H319	Causes serious eye irritation.
	H336	May cause drowsiness or dizziness.
	Precautionary statements	
	P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
	P233	Keep container tightly closed.
	P280	Wear eye protection.
	P337 + P313	If eye irritation persists: Seek medical advice/attention.
	P370 + P378	In case of fire: Use powder/sand/carbon dioxide extinguisher to extinguish.
	P403 + P235	Store in a well-ventilated area. Keep cool.

MT-002 — Hydrolysis Reagent

Parameter	Hazard statements	
  	H314	Causes severe skin burns and eye damage.
	H370	Causes damage to organs.
	H301 + H311 + H331	Toxic if swallowed, in contact with skin, or if inhaled.
	Precautionary statements	
	P260	Do not breathe dust/fume/gas/mist/vapours/spray.
	P280	Wear protective gloves/protective clothing/eye protection/face protection.
	P301 + P330 + P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
	P303 + P361 + P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
	P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
	P310	Immediately call a doctor.

MT-003 — Neutralization Reagent

Parameter	Hazard statements	
	H315	Causes skin irritation.
	H319	Causes serious eye irritation.
	H335	May cause respiratory irritation.
	Precautionary statements	
	P261	Avoid breathing mist/vapours/spray.
	P264	Wash hands and exposed parts of the body thoroughly after handling.
	P280	Wear protective gloves.
	P312	Call a POISON CENTRE if you feel unwell.
	P321	Specific treatment (see supplemental first aid instructions on this label).
	P337 + P313	If eye irritation persists: Get medical advice/attention.

MT-004 — Extraction Reagent

Parameter	Hazard statements	
 	H302	Harmful if swallowed.
	H315	Causes skin irritation.
	H319	Causes serious eye irritation.
	H331	Toxic if inhaled.
	H351	Suspected of causing cancer.
	H361d	Suspected of damaging the unborn child.
	H372	Causes damage to organs through prolonged or repeated exposure.
	Precautionary statements	
	P260	Do not breathe vapours.
	P264	Wash hands and exposed parts of the body thoroughly after handling.
	P280	Wear protective gloves.
	P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
	P311	Call a doctor.
	P403 + P233	Store in a well-ventilated place. Keep container tightly closed.

MT-005 — Dilution Reagent

Parameter	Hazard statements	
  	H225	Highly flammable liquid and vapour.
	H370	Causes damage to organs.
	H301 + H311 + H331	Toxic if swallowed, in contact with skin, or if inhaled.
	Precautionary statements	
	P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
	P280	Wear protective gloves/protective clothing.
	P301+P310	IF SWALLOWED: Immediately call a doctor.
	P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
	P308+P311	IF exposed or concerned: Call a doctor.
	P311	Call a doctor.
	P403 + P235	Store in a well-ventilated place. Keep cool.

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Website: www.dimedical.pl

Additional information

The manufacturer holds a certified quality management system compliant with **ISO 13485**.

Device identification

Element	Content
Device name	M-Typer Reagent Kit (Set 1 + Set 2)
Kit catalog number	MT-KIT-001
Set 1 catalog number	MT-KIT-SET1
Set 2 catalog number	MT-KIT-SET2
Batch / Lot number	Stated on the packaging label
Expiry date	Stated on the packaging label
IFU version number	IFU-M-TYPER-001-v1.0
IFU issue date	[to be completed upon approval]
Device classification	In vitro diagnostic medical device (IVD)
Device class	IVDR Class C

Section B: Reporting of Serious Incidents

A serious incident is an event related to the use of a medical device that has led or could lead to:

- death of the patient, user, or another person,
- serious deterioration of the health of the patient, user, or another person,
- a serious threat to public health.

In the event of a confirmed or suspected serious incident related to the ***M-Typer Reagent Kit***, the event must be reported immediately to the manufacturer and to the relevant national competent authority.

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In the event of external corrective safety actions (**FSCA**) taken by the manufacturer, users of the device will be informed by means of a Field Safety Notice (**FSN**).

Section C: Explanation of Symbols

The following symbols may appear on the packaging, label, or in these Instructions for Use of the *M-Typer Reagent Kit*.

Symbol	Meaning
	Manufacturer
	Country of manufacture
	Date of manufacture
	Use by (expiry date)
	Batch / Lot number
	Catalog number
	In vitro diagnostic medical device
	Contents sufficient for 100 tests
	Temperature limitation for storage
	Protect from light
	Contains hazardous substances (carcinogenic, mutagenic, or toxic to reproduction — CMR)
	Do not use if packaging is damaged, and consult the instructions for use
	Caution
	Consult instructions for use
	Unique Device Identifier (UDI)
	CE marking of conformity